Saccharoflex 2000

Electronic whiteness-degree meter

- Degree of whiteness determination of crystal sugar
- Value according to ICUMSA method GS2-11 and GS2-13
- Easy calibration by Brunswick Color Standards
- Designed for permanent use

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Saccharoflex 2000

The SCHMIDT+HAENSCH electronic reflectance-meter for the automatic color grading of crystallized sugar according to the recommendations of ICUMSA, method GS2-11 (1994) and GS2-13 (1994).

With this instrument the reflectance factors at 495 nm and 620 nm are measured and the ratio is indicated in terms of color type numbers. The measurement is done by use of Taylor's law and an integrated sphere. To decrease measurement instabilities caused by inhomogeneity of the sample (e.g. different grain size) a great sample surface was chosen.

Only two standards of the Brunswick color type, preferably 0 + 6 are necessary to calibrate the Saccharoflex 2000.

An adjustment can be carried out by the keyboard. The Saccharoflex 2000 allows measurements up to 19.99 indicating units without further adjustments. LCD display with 4-button keyboard makes handling easy and standard serial interface RS 232 C ensures proper printing.

The adjustable working voltage provides a high service life of the lighting unit and by this large replacement intervals. In addition the electric bulb is run by a power saving switch which prevents the Saccharoflex 2000 from heating up, so that a continuous operation is possible without any problems.

The consistent use of high integrated circuits, the little heating and the easily detachable lighting unit ensure that the Saccharoflex 2000 is a reliable and easy maintainable measurement system which corresponds to the high SCHMIDT+HAENSCH quality standard.

Technical data:

Keyboard:	4-button keyboard
Display:	LCD, one line, 30 character
Measuring range:	0 to 19.99 color units
Resolution:	0.01 CTU and 0.05 CTU
Wavelength:	495 + 620 nm
Calibration:	Standards 0 + 6 (Brunswick)
Data output:	Serial interface RS 232 C
Sample compartment:	Round table (standard ICUMSA)
Dimensions, weight:	22 x 26 x 48 cm, 9 kg
Accessories (necessary):	Color standard 0 + 6, each 200 g

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SCHMIDT + HAENSCH

Opto-electronic measuring device since 1864

ISO 9001:2000

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SCHMIDT+HAENSCH GmbH & Co.